

# VITA-D-CHLOR™

Use the best dechlorination chemical with the best dechlorination devices.

## DECHLOR DEMON™

### Performance Data

#### Chlorine Levels

0–50 ppm

#### Flush Rate

150–900 gpm

#### Setup Time

2–7 mins.



### Procedure (Treatment up to 3 ppm chlorine)

1. Setup the Dechlor Demon to Hydrant (Hose Monster diffuser should be used for safety)
2. Add Vita-D-Chlor tablets to mixing tank
3. Fill tank with water and connect to Dechlor Demon
4. Open outlet ball-valve fully and inlet ball-valve halfway
5. Begin flush and adjust ball valves for greatest efficiency and test for complete dechlorination

### Benefits/Features

- Dechlorinate while flow testing or flushing
- Minimal friction loss
- Does not stress hydrant nozzles
- Available in 2.5", 4" and 4.5" models
- Durable construction
- Small, compact design



Manufactured by Hydro Flow Products (Hose Monster Co.)

### VITA-D-CHLOR™

- NSF Granular
- One gram neutralizes 100 gallons of water with 1 ppm chlorine
- Soluble
- Safe for operators & the environment

#### AVAILABLE UNIT SIZES:

- 2.5 kg
- 12 kg
- 25 kg resealable poly pails

### VITA-D-CHLOR™ Neutral

- NSF Granular
- pH Neutral
- Very Soluble
- Safe for operators & the environment

#### AVAILABLE UNIT SIZES:

- 2.5 kg
- 12 kg
- 25 kg resealable poly pails

### VITA-D-CHLOR™ Tablets

- 2.5 inch tablet
- Easy to use
- 100% organic and water soluble
- Safe for operators & the environment

#### AVAILABLE UNIT SIZES:

- 20 tablets
- 40 tablets
- 35 lb (approx. 140 tablets)



**INTEGRA**  
chemical

■ Phone 800.322.6646 ■ Website [www.Vita-D-Chlor.com](http://www.Vita-D-Chlor.com)